

QuickGuide: RealFast™ Variant Detection on LightCycler® 480II

Setup for Variant Detection RealFast™ Assays:

- Open the **LightCycler® 480 software** (QuickGuide is based on version 1.5.1) and login with your username and password.
 - Choose **New Experiment**.
- Define your PCR program in the **Run Protocol** tab.
 - Select **Dual Color Hydrolysis Probe / UPL Probe** as Detection Format.
 - Select a **reaction volume** of 20 µl and setup the **Program**:
 - Add programs in the **Program Name** window with "+" and edit **Cycles** and **Analysis Mode**. To edit **Acquisition Mode** and **Hold**, click on the corresponding step in the **Program Name** window and change parameters or add steps ("+") in the **Temperature Targets** window.
 - Setup the PCR program according to the Instructions for Use of your RealFast™ Assay.
- Click on **Subset Editor** button on the left side of the window.
 - Press "+" to create a new subset and rename your subset.
 - Select wells in the grid and press **Apply**.
- Click on **Sample Editor** button on the left side of the window.
 - **Select Workflow > Abs Quant.**
 - **Select Filter Combinations > 456-510 nm (FAM) and 533-580 nm (HEX).**
 - Choose your **Subset** of Samples.



- Define your **No Template Control (NTC)**:

- **Select Samples** field: select well by mouse-click or two wells by ctrl + mouse click.
- **Edit Abs Quant Properties** field: Type **NTC** in the **Sample Name** field and press **Enter**. Choose **Negative Control** as **Sample Type**.

| Pos | Filter combination | Color | Repl Of | Sample Name | Quantification Sample Type | Concentration |
|------|--------------------|-------------------------------------|---------|-------------|----------------------------|---------------|
| ▶ A1 | FAM (465-510) | ■ | A1 | NTC | Negative C | |
| A1 | VIC / HEX / | ■ | A1 | NTC | Negative Con | |
| A2 | FAM (465-510) | ■ | A1 | NTC | Negative Con | |
| A2 | VIC / HEX / | ■ | A1 | NTC | Negative Con | |

- Define your **Samples**:

- **Select Samples** field: select a well by mouse click.
- **Edit Abs Quant Properties** field: Type the name of your first sample in the corresponding field and press **Enter**. Check **Unknown**.
- Define the rest of your samples.

| Pos | Filter combination | Color | Repl Of | Sample Name | Quantification Sample Type | Concentration |
|------|--------------------|--|---------|-------------|----------------------------|---------------|
| ▶ A3 | FAM (465-510) | ■ | A3 | 3 | Unknown | |
| A3 | VIC / HEX / | ■ | A3 | 3 | Unknown | |
| A4 | FAM (465-510) | ■ | A3 | 3 | Unknown | |
| A4 | VIC / HEX / | ■ | A3 | 3 | Unknown | |

- **Save** or **export** your experiment by pressing the corresponding button:



- Load your samples and start the experiment.




Analysis of Variant Detection RealFast™ Assays using “Abs Quant/Fit Points”:

- Open the **LightCycler® 480 software** and login with your username and password. The **Overview** window appears.
- Click on  or choose **Navigator** in the flip-window on the top left.
- Choose an experiment from the data bank or
- import an experiment located outside the data bank by pressing **Import**.
- After the file is loaded the **Summary** window of your experiment is displayed.
- Press the **Analysis** button to reach the analysis window.

- Within the **Create New Analysis** field choose **Abs Quant/Fit Points**.
- A pop-up window will be launched. If applicable select a **Subset** of samples or analyze **All Samples** in case your plate contains only one type of assay. Give a **Name** to your analysis. Press the **OK** button

- *Optional:* Press **Color Compensation** and choose **In Use** for color compensation of FAM (510) and VIC (580). Within the pop-up window select a color compensation and press the **OK** button.
- Press **Filter Comb 465-510**, choose **FAM (465-510)** in the pop-up window and press the **OK** button.
- Within the **Noise Band** field adjust the threshold manually. Set the threshold in the exponential phase of the amplification curves above the background.

| Abs Quant results | | | | |
|---|--------------------------------------|-----|----------|--------|
| <input type="checkbox"/> Positive <input type="checkbox"/> Negative <input type="checkbox"/> Standard | | | | |
| Samples | | | | |
| Include | Color | Pos | Name | Cp |
| <input checked="" type="checkbox"/> | ■ | A4 | B*070201 | |
| <input checked="" type="checkbox"/> | ■ | B4 | B*070201 | |
| <input checked="" type="checkbox"/> | ■ | C4 | B*2701 | 24, 45 |
| <input checked="" type="checkbox"/> | ■ | D4 | B*2701 | 24, 59 |
| <input checked="" type="checkbox"/> | ■ | E4 | H2O | |
| <input checked="" type="checkbox"/> | ■ | F4 | H2O | |



- Press the **Calculate** button. Choose **Abs Quant Results** in the chart bottom-left and review your samples. Samples positive for the targeted variant show an amplification in the FAM channel and will be color-coded.
- To review proper amplification of all samples, press **Filter Comb 465-510**, choose **VIC / HEX / Yellow 555 (533/580)** in the pop-up window and press the **OK** button.
- Press the **Calculate** button and review your results in the chart bottom-left. All genomic DNA samples and Controls show an amplification in the HEX channel and are therefore marked as positive.


» **Note:** Set threshold value for the FAM channel just above the background fluorescent signal generated by the e.g. HLA-B27 Negative Control. «

- **Save** or **export** your data by pressing the corresponding button.  or 
- After saving your data you can customize and generate a report via the **Report** button.

Optional:

Analysis of Variant Detection RealFast™ Assays using Abs Quant/2nd Derivative Max:


- Open the **LightCycler® 480 software** and login with your username and password. The **Overview** window appears.
- Click on  or choose **Navigator** in the flip-window on the top left.
 - Choose an experiment from the data bank or
 - import an experiment located outside the data bank by pressing **Import**.
- After the file is loaded the **Summary** window of your experiment is displayed.
- Press the **Analysis** button to reach the analysis window.
 - Within the **Create New Analysis** field choose **Abs Quant/2nd Derivative Max**.
 - A pop-up window will be launched. If applicable select a **Subset** of samples or analyze **All Samples** in case your plate contains only one type of assay. Give a **Name** to your analysis. Press the **OK** button.




Create new analysis

| | |
|---------------|--|
| Analysis Type | * Abs Quant/2nd Derivative Max |
| Subset | * HLA B57 |
| Program | * amplification |
| Name | * Abs Quant/2nd Derivative Max for HLA B57 |

- *Optional:* Press **Color Compensation** and choose **In Use** for color compensation of FAM (510) and VIC (580). Within the pop-up window select a color compensation and press the **OK** button.
- Press **Filter Comb 465-510**, choose **FAM (465-510)** in the pop-up window and press the **OK** button.
- Press the **Calculate** button. Choose **Abs Quant Results** in the chart bottom-left and review your samples. Samples positive for the targeted variant show an amplification in the FAM channel and will be color-coded.



| Include | Co... | Pos | Name | Cp |
|-------------------------------------|-------|-------|--------|-------|
| <input checked="" type="checkbox"/> | | Green | A1 NTC | |
| <input checked="" type="checkbox"/> | | Green | A2 NTC | |
| <input checked="" type="checkbox"/> | | Green | A5 4 | |
| <input checked="" type="checkbox"/> | | Green | A6 4 | |
| <input checked="" type="checkbox"/> | | Red | A9 9 | 25,89 |
| <input checked="" type="checkbox"/> | | Red | A10 9 | 25,81 |



- To review proper amplification of all samples, press **Filter Comb 465-510**, choose **VIC / HEX / Yellow 555 (533/580)** in the pop-up window and press the **OK** button.
 - Press the **Calculate** button and review your results in the chart bottom-left. All DNA samples and the Positive Control show an amplification in the HEX channel and are therefore marked as positive.
- » **Note:** Set threshold value for the FAM channel just above the background fluorescent signal generated by the e.g. HLA-B27 Negative Control. «

- **Save** or **export** your data by pressing the corresponding button.  or 
- After saving your data you can customize and generate a report via the **Report** button.