


QuickGuide: RealFast™ Genotyping on Rotor-Gene® 6000

Setup for Genotyping Assays:

- Launch the **Rotor-Gene® 6000 software** (QuickGuide is based on version 1.7).
- From the upper menu press **New**.
 - Within the New Run window select the **Advanced** tab.
 - Choose **Two Step** and press **New**.
- The **New Run Wizard** window will appear.
 - Select the **Rotor Type** and tick the box **Locking Ring Attached**. Press **Next**.
 - Define a Reaction Volume of **20 µl**. *Optional: enter an operator name and run specific notes.* Press **Next**.
 - Click **Edit Profile**, enter the following PCR program and confirm with **OK**:

For **singleplex** RealFast™ Assays use:

	Cycles	Temperature [C°]	Duration [mm:ss]	Acquiring
Hold	-	95	03:00	-
Cycling	40	95	00:15	-
		36-well rotor: 56 72-well rotor: 60	01:00	Acquiring to Cycling A on Green and Yellow




For **multiplex** RealFast™ Assays use:

	Cycles	Temperature [C°]	Duration [mm:ss]	Acquiring
Hold	-	95	03:00	-
Cycling	40	95	00:15	-
		36- well & 72-well rotor: 60	01:00	Acquiring to Cycling A on Green, Yellow, Red, Orange

- Click **Gain Optimisation**
- Select **Green** from the drop down menu **Channel Settings**. Press **Add**
Keep the default values in **Auto-Gain Optimisation Channel Settings** and confirm.

Select **Yellow** from the drop down menu **Channel Settings**. Press **Add**
Keep the default values in **Auto-Gain Optimisation Channel Settings** and confirm.



Repeat this procedure with **Orange** and **Red** for **multiplex RealFast™ Assays**.
 Tick the box **Perform Optimisation Before 1st Acquisition**. Close window.

» **Note:** *In case you carry out different RealFast™ Assays within one run, it is advisable to do gain optimisation individually for each of the assays. Create new channels, e.g. FAM1, HEX1,.....for data acquisition.* «

Setting	Value
FAM1 Gain	5
Green Gain	5
HEX1 Gain	5
Yellow Gain	5
Auto-Gain Optimisation	Before First Acquisition
Rotor	72-Well Rotor
Sample Layout	1, 2, 3, ...
Reaction Volume (in microliters)	20

➤ Proceed with **Next**.

- Press **Start Run**.
- Define a file name and press **Save**.
- The **New Run Wizard** window will pop up, where the sample **Names** and **Types** (NTC, Unknown and Positive Controls) can be defined.
- Press **Finish** and the **Run Progress** window will appear.

Analysis of Genotyping Assays using Allelic Discrimination:

- Launch the **Rotor-Gene® 6000 software**.
- In the upper menu press **Open** and load your genotyping data file.
- The software displays the **Raw Channel (Cycling A. Green)** window.
- Press the **Analysis** button from the upper menu and choose **Allelic Discrimination**.
 ➤ Select **Cycling A. Green** and **Cycling A. Yellow** and press **Show**.

- In the **Allelic Discrimination Analysis** window press **Genotypes**.


- According to your assay, select the appropriate channels for **wild type (WT)** (mostly **Cycling A. Yellow**), **mutant (MUT)** (mostly **Cycling A. Green**) and **Heterozygous** (**Cycling A. Green** and **Cycling A. Yellow**). Confirm with **OK**.

Genotype	Reacting Channels	
Wild Type		Cycling A.Yellow
Heterozygous	Cycling A.Green	Cycling A.Yellow
Mutant	Cycling A.Green	

- Press **Dynamic Tube**.

*Optional: If necessary, use either **Slope Correct** or **Ignore First** for optimization of amplification curves, whatever fits best.*

- Set the **Discrimination Threshold** as following:
 ➤ Click on the threshold button  in the lower right panel. In the analysis graph click and drag the red threshold line above the background signals of the positive controls.
- Results are shown in the **Allelic Discrimination Results** window.
- To generate a report press **Reports** from the upper menu and select the relevant report in the **Report Browser** to be shown and printed.



» **Note:** *With **multiplex RealFast™ Assays** the two markers (corresponding to HEX/FAM and Cy5/ROX, respectively) have to be analyzed **one after the other** !* «

- Repeat the analysis for the other marker by pressing **Analysis** → **Allelic Discrimination**. Then select **Cycling A. Orange** and **Cycling A. Red** and press **Show**.

- Proceed as described above. For most assays **wild type (WT)** corresponds to **Cycling A. Red** and **mutant (MUT)** to **Cycling A. Orange**.

Genotype	Reacting Channels	
Wild Type		Cycling A.Red
Heterozygous	Cycling A.Orange	Cycling A.Red
Mutant	Cycling A.Orange	

» **Note:** *For assignment of fluorophores to hydrolysis probes see **Instructions for Use of the respective RealFast™ Assay**.* «

Analysis of Genotyping Assays using Scatter Plot:

- Launch the **Rotor-Gene® 6000 software**.
- In the upper menu press **Open** and load your genotyping data file.
- The software displays the **Raw Channel (Cycling A. Green)** window.
- Press the **Analysis** button from the upper menu and choose **Scatter Graph Analysis**.
 - Select **Cycling A. Green** and **Cycling A. Yellow** and press **Show**.


- In the **Scatter Analysis** window press **Genotypes**.
 - According to your assay, select the appropriate channels for **wild type (WT)** (mostly **Cycling A. Yellow**), **mutant (MUT)** (mostly **Cycling A. Green**) and **heterozygous (Cycling A. Green and Cycling A. Yellow)**. Confirm with **OK**.

Genotype	Reacting Channels	
Wild Type		Cycling A.Yellow
Heterozygous	Cycling A.Green	Cycling A.Yellow
Mutant	Cycling A.Green	

- Press **Dynamic Tube**.

*Optional: If necessary, use either **Slope Correct** or **Ignore First** for optimization of amplification curves, whatever fits best.*

- Data points are shown in the **Scatter Analysis Graph** window.
 - Draw a rectangle around the data points in each quarter > define the corresponding genotypes as **Wild Type**, **Mutant** and **Heterozygous**, respectively, and the No-Template Control as **None**.



- Samples and their corresponding genotypes are shown in the **Scatter Analysis Result** window.
- To generate a report press **Reports** from the upper menu and select the relevant report in the **Report Browser** to be shown and printed.

» **Note:** With **multiplex RealFast™ Assays** the two markers (corresponding to **HEX/FAM** and **Cy5/ROX**, respectively) have to be analyzed **one after the other** ! «

- Repeat the analysis for the other marker by pressing **Analysis** → **Allelic Discrimination**. Then select **Cycling A. Orange** and **Cycling A. Red** and press **Show**.

- Proceed as described above. For most assays Wild Type corresponds to **Cycling A. Red** and Mutant to **Cycling A. Orange**.

Genotype	Reacting Channels	
Wild Type		Cycling A.Red
Heterozygous	Cycling A.Orange	Cycling A.Red
Mutant	Cycling A.Orange	

» **Note:** For assignment of fluorophores to hydrolysis probes see *Instructions for Use of the respective RealFast™ Assay*. «