


## QuickGuide: RealFast™ CNV on StepOne™

### Setup of Relative Quantitation Assays:

- Open the StepOne™ Software (QuickGuide is based on version 2.3) and click **New Experiment** and **Advanced Setup**.
  - In **Setup > Experiment Properties** select:
    - Instrument: **StepOne™** (48 Wells)
    - Type of experiment: **Quantitation – Comparative C<sub>T</sub> (ΔΔC<sub>T</sub>)**
    - Reagents: **TaqMan® Reagents**
    - Ramp speed: **Standard**
  - In **Setup > Plate Setup** go to **Define Targets and Samples**:
    - **Define Targets** in the corresponding field:
      - Provide a name for your gene of interest and choose **FAM** as **Reporter** and **NFQ-MGB** as **Quencher**.
      - Add a **New Target** by pushing the corresponding button.
      - Type **EC** (endogenous control) as **Target Name** and choose **VIC** as **Reporter** and **NFQ-MGB** as **Quencher**.
- | Target Name             | Reporter | Quencher | Color |
|-------------------------|----------|----------|-------|
| Gene of interest        | FAM      | NFQ-MGB  | Red   |
| EC (endogenous control) | VIC      | NFQ-MGB  | Blue  |
- **Define Samples** in the corresponding field:
      - Type **Calibrator** in the field for **Sample Names**. This represents the positive control which is included in the assay kit.
      - **Add New Sample(s)** by pushing the corresponding button and rename the field(s) according to the sample(s) you want to analyze.
  - In **Setup > Plate Setup** go to **Assign Targets and Samples**
    - Define the **Negative Control Template**:
      - Select a replicate of three wells by ctrl-click.
      - Within the field **Assign Target(s) to the selected wells** check boxes for the gene of interest (e.g. CYP21A2) and **EC**. Click on the button **N** (Negative Control) in **Task**.
    - Define your **Calibrator**:
      - Select a replicate of three wells by ctrl-click.
      - Within the field **Assign Target(s) to the selected wells** check boxes for the gene of interest (e.g. CYP21A2) and **EC**. Click on the button **U** (Unknown) in **Task**.
      - Check the box for the **Calibrator** within the field **Assign Sample(s) to the selected wells**.
    - Within the field **Select relative quantitation settings** choose **Calibrator** as your **Reference Sample** and **EC** as **Endogenous Control**.
    - Within the field called **Select the dye to use as the passive reference** select **ROX**.
    - Define your **Samples**:
      - Select a replicate of three wells by ctrl-click.
      - Within the field **Assign Target(s) to the selected wells** check boxes for the gene of interest (e.g. CYP21A2) and **EC**. Click on the button **U** (Unknown) in **Task**.
      - Check the box for the **Sample** you wish to assign within the field called **Assign Sample(s) to the selected wells**.
  - In **Setup > Run Method** go to **Graphical View**
    - Select a reaction volume of **20 µl**
    - Define your PCR program:
      - *optional: include Pre-PCR Read*
      - Holding Stage: **10 min at 95°C**
      - Cycling Stage: **40 cycles 15 sec at 95°C and 1 min at 60°C**. Make sure **Data Collection On** is enabled
  - Load your reaction plate into the StepOne™ instrument and press **START RUN** (green button).

Assign	Target	Task
<input checked="" type="checkbox"/>	Gene of interest	<input type="checkbox"/> <input type="checkbox"/>
<input checked="" type="checkbox"/>	EC (endogenous ...)	<input type="checkbox"/> <input type="checkbox"/>


Mixed  Unknown  Negative Control



## Analysis of Relative Quantitation Assays:

After completing a run or after opening a genotyping data file the software displays the Experiment Menu **Analysis:**

- Results automatically appear in the **Amplification Plot**.
  - Adjust the **Plot Settings** to  $\Delta Rn$  vs Cycle (Plot Type), **Linear** or **Log** (Graph Type), **Target** (Color)




Plot Settings

Plot Type:  $\Delta Rn$  vs Cycle    Graph Type: Linear    Color: Target

Save current settings as the default

- Tick the box for **Show Threshold** in the **Options** field.
- Press the **Analysis Settings** Button and go to **C<sub>T</sub> Settings**.
  - Adjust the **Threshold** according to the settings in the **Assay Description** and press the button Apply Analysis Settings.
- Select individual replicates in the **View Plate Layout** field and review your samples.
  - The interval between the curve for the gene of interest and for the **Endogenous Control (EC)** is related to the copy number variation.



Amplification Plot of the Calibrator sample. Linear (left) and log (right) graph type.

Example:

- Go to **View Well Table**.
  - Press the **Show in Table** button and customize the table.
  - Review the **Relative Quantities (RQ)** and define the CNV status of your samples according to the Assay Description.
  - Go to **Gene Expression** (left) and select **RQ vs Sample** in the **Plot settings**. The relative quantities of each sample are displayed as bar chart.
- To print a report click **Print Report** in the upper menu bar:
  - Select data for the report according to your needs.