StripAssay[®] Detection Reagents 48

REF

CS-017

IVD

48 Tests

2-8°C

CE

1. DNAT (<i>blue cap</i>)	1.5 ml	⚠ Warning
2. Hybridization Buffer (<i>white cap</i>)	65 ml	
3. Wash Solution A (<i>white cap</i>)	200 ml	
4. Conjugate Solution	65 ml	
5. Wash Solution B	200 ml	
6. Color Developer	65 ml	⚠ Warning
7. Instructions For Use	1	

ViennaLab Diagnostics GmbH

Gaudenzdorfer Guertel 43-45

A-1120 Vienna, Austria

Phone: (+43-1) 8120156-0

info@viennalab.com

www.viennalab.com

Instructions for use

I. INTENDED USE

StripAssay® Detection Reagents 48 are to be used in combination with any ViennaLab StripAssay®. These are auxiliary reagents for the hybridization on automated instruments for up to 48 Teststrips. For a detailed assay procedure and interpretation of results follow the instructions for use of the respective StripAssay®. *For human in vitro diagnostics.*

III. KIT COMPONENTS

See list of all kit components on page I.

Warning: DNAT contains 1.6% NaOH.

H315: Causes skin irritation

H319: Causes serious eye irritation

P280: Wear protective gloves/protective clothing/eye protection/face protection

P337 + P313: If eye irritation persists: Get medical advice/attention

Warning: Color Developer contains ≤0.4% maleic acid.

H317: May cause an allergic skin reaction

P280: Wear protective gloves/protective clothing/eye protection/face protection

P302 + P352: If on skin: wash with plenty of water

P333 + P313: If skin irritation or rash occurs: get medical advice /attention

Amplification Mix, Taq Dilution Buffer, Conjugate Solution, Wash Solution B contain 0.05% NaN₃. Conjugate Solution contains streptavidin-alkaline phosphatase. Color Developer contains nitro blue tetrazolium (NBT) and 5-bromo-4-chloro-3-indolyl phosphate (BCIP).

Store all reagents at 2-8°C when not in use !

VII. QUALITY CONSIDERATIONS

- Use of the StripAssay for human *in vitro* diagnostics needs to be limited to appropriately trained personnel.
- Do not use StripAssay components beyond the expiration date printed on the outside of the kit box. Do not mix reagents from different lots.
- Avoid microbial contamination and cross-contamination of reagents or samples by using sterile disposable pipette tips throughout. Do not interchange bottle caps.

VIII. SAFETY

- Do not drink, eat, smoke, or apply cosmetics in designated work areas. Wear laboratory coats and disposable gloves when handling specimens and kit reagents. Wash hands thoroughly afterwards.
- Handle specimens as if capable of transmitting infectious agents. Thoroughly clean and disinfect all materials and surfaces that have been in contact with specimens. Discard all waste associated with clinical specimens in a biohazard waste container.
- Avoid contact of DNAT with skin, eyes, or mucous membranes. If contact does occur, immediately wash with large amounts of water. If spilled, dilute with water before wiping dry.
- Adhere to all local and federal safety and environmental regulations which may apply.

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