MTHFR mpx RealFast™ Assay

ViennaLab Diagnostics GmbH

Gaudenzdorfer Guertel 43-45

A-1120 Vienna, Austria

Phone: (+43-1) 8120156-0

info@viennalab.com

www.viennalab.com

REF 7-165 / 7-168 Σ 100 / 32 reactions

-30°C / -15°C

1. Intended Use

The MTHFR mpx RealFast™ Assay is a fast and accurate multiplex real-time PCR test for the simultaneous detection of the 677C>T and 1298A>C mutations in the *methylenetetrahydrofolate reductase (MTHFR)* gene. The mutations are associated with decreased enzyme activity, which leads to hyperhomocysteinemia and to toxic side-effects of methotrexate therapy. The kit is designed to identify patients suspected to have an increased risk for cardiovascular diseases or intolerance to methotrexate. The qualitative assay discriminates the three possible genotypes for each allele in a human DNA extract.

Reference sequence: NG_013351.1 g.14783C>T and g.16685A>C; dbSNP: rs1801133 and rs1801131.

2. Introduction

Elevated levels of homocysteine, an intermediary product of the methionine metabolism, are an established risk factor for atherosclerosis and arterial thrombosis. Hyperhomocysteinemia ultimately leads to endothelial dysfunction with associated platelet activation and thrombus formation. Mild to moderate forms can be caused by homozygosity for a common 677C>T point mutation in the coding region of the *MTHFR* gene. The second common MTHFR mutation 1298A>C also contributes to reduced enzyme activity, especially when occurring simultaneously with 677C>T. Fasting homocysteine levels are significantly higher in individuals heterozygous for both substitutions compared with individuals carrying only the heterozygous 677C>T mutation. Hyperhomocysteinemia may affect methotrexate sensitivity and contribute to toxicity, as MTHFR is an important enzyme in maintaining cellular folate pools. Thus, homozygosity for 677C>T or compound heterozygosity for 677C>T/1298A>C conveys a significantly higher risk for negative side-effects of methotrexate medication.

3. Kit Contents

100 / 32 Rxn

RealFast™ 2x mpx Probe Mix	1 vial white cap	1000 / 320 µl	The RealFast™ 2x mpx Probe Mix comprises HotStart Taq DNA polymerase and dNTPs in an optimized buffer system. The MTHFR mpx Assay Mix consists of <i>MTHFR</i> gene-specific primers and four allele-specific, dual-labeled hydrolysis probes. Controls representing wild type (WT-Control) and homozygous mutant (MUT-Control) genotypes are supplied with the kit.
MTHFR mpx Assay Mix	1 vial purple cap	550 / 550 µl	
MTHFR mpx WT-Control	1 vial green cap	75 / 75 µl	
MTHFR mpx MUT-Control	1 vial red cap	75 / 75 µl	

The kit contains reagents for 100 / 32 reactions in a final volume of 20 µl each.

4. Storage and Stability

The MTHFR mpx RealFast™ Assay is shipped on cooling blocks. On arrival, store the kit at -30 to -15°C. Alternatively, store at 2 to 8°C for short-term use within one month. The kit withstands up to 20 freeze/thaw cycles with no loss of activity. Avoid prolonged exposure to intense light. If stored correctly, the kit will retain full activity until the expiration date indicated on the label.

5. Product Description

5.1. Principle of the Test

The test is based on the fluorogenic 5' nuclease assay, also known as TaqMan® assay. Each reaction contains two gene-specific primer pairs which amplify a 144 bp and 185 bp fragment of the *MTHFR* gene, as well as four dual-labeled, allele-specific hydrolysis probes which specifically hybridize to the target sequences of the amplified fragments. The proximity of the 5'-fluorescent reporter and 3'-quencher dye on intact probes prevents the reporter from fluorescing. During the extension phase of PCR the 5' – 3' exonuclease activity of the Taq DNA polymerase cleaves the 5'-fluorescent reporter from the hybridized probe. The physical separation of the fluorophore from the quencher dye generates a fluorescent signal in real-time, which is proportional to the accumulated PCR product.

MTHFR probe	Fluorophore	Channel
677C>T mutant	FAM	520 nm
677C>T wild type	HEX	556 nm
1298A>C mutant	ROX	605 nm
1298A>C wild type	Cy5	670 nm

In normal samples the **wild type probes** generate a strong fluorescence signal in the HEX or Cy5 channel and no or only a baseline signal in the FAM or ROX channel. Vice versa, in homozygous mutant samples the hybridized **mutant probes** generate a strong fluorescence signal in the FAM or ROX channel and no or only a baseline signal in the HEX or Cy5 channel. In heterozygous samples, both wild type and mutant probes bind to the amplicons and generate intermediate signals in the respective channels.

5.2. Real-time PCR Instrument Compatibility

The MTHFR mpx RealFast™ Assay is validated for use with the AB 7500 Fast instrument.

The kit is compatible with various common real-time PCR instruments capable of recording FAM, HEX, Cy5 and ROX fluorescence:

- ✓ AB 7500 Fast (Applied Biosystems®)
- ✓ CFX96™ (Bio-Rad)
- ✓ LightCycler® 480 (Roche)
- ✓ MIC qPCR Cycloer (bms)
- ✓ Rotor-Gene® 6000 (Qiagen)

» **Note:** RealFast™ Genotyping QuickGuides for setting up and analyzing experiments on different types of instruments can be downloaded from www.viennalab.com.

The kit is **not suitable** for use with real-time PCR instruments requiring ROX for normalization of data (e.g. Applied Biosystems® instruments: StepOne™, 7300, 7900/7900HT) or for instruments without appropriate fluorescence detection channels.

5.3. Assay Performance Specifications

Determination of **sensitivity** was performed on 69 positive 677C>T alleles and 69 positive 1298A>C alleles, both tested with a CE-marked reference kit. The MTHFR mpx RealFast™ Assay correctly determined all positive MTHFR alleles, which equaled a true positive rate of 100%.

Determination of **specificity** was performed on 113 negative 677C>T alleles and 113 negative 1298A>C alleles, both tested with a CE-marked reference kit. The MTHFR mpx RealFast™ Assay correctly determined all negative MTHFR alleles, which equaled a true negative rate of 100%.

Limit of detection: 0.2 ng genomic DNA (per reaction). Recommended DNA concentration: 2 to 20 ng/µl genomic DNA.

6. Materials Required but not Supplied

Real-time PCR instrument with FAM (520 nm), HEX (556 nm), ROX (605 nm) and Cy5 (670 nm) filters, instrument-compatible reaction vessels, disposable powder-free gloves, vortexer, mini-centrifuge for 2.0 ml tubes, tube racks, set of calibrated micropipettes (0.5 – 1000 µl), sterile tips with aerosol-barrier filter, molecular grade water, DNA extraction system, freezer, biohazard waste container.

7. Experimental Protocol

7.1. DNA Extraction

DNA extraction reagents are **not supplied** with the kit.

DNA isolated from various specimens (e.g. whole peripheral blood, dried blood spots, buccal swabs or saliva) can be used. Ensure extracted DNA is suitable for amplification in terms of concentration, purity and integrity.

For accurate genotype calling, the DNA amount per reaction should be within the range of 10 to 100 ng for all samples.

7.2. PCR Controls

Always include a **No Template Control (NTC)** in each experiment to confirm absence of potential contamination. It is advisable to run the NTC (use PCR-grade water instead of DNA) in duplicate.

Always include the MTHFR mpx **WT-Control** and MTHFR mpx **MUT-Control** as positive reference signals for your unknown samples. Some real-time PCR software, e.g. AB 7500 Fast, requires signals for all three possible genotypes for correct allelic discrimination. In order to obtain a heterozygous control (HET-Control), mix an aliquot of WT-Control and MUT-Control in a ratio of 1:1.

» **Note:** *WT- and MUT-Controls are potential sources of contamination. Make sure to handle them carefully.* «

7.3. Preparation of MTHFR mpx RealFast™ Master Mix:

Gently vortex and briefly centrifuge all solutions after thawing. Set up PCR at room temperature. Prepare sufficient **Master Mix** for all your reactions (N samples + positive controls + negative controls) plus at least one additional reaction to compensate for pipetting inaccuracies:

Component	per reaction	e.g. 24+1 reactions
RealFast™ 2x mpx Probe	10 µl	250 µl
MTHFR mpx Assay Mix	5 µl	125 µl
Master Mix	15 µl	375 µl

Dispense **15 µl Master Mix** into each well. Add **5 µl purified DNA** or **Control** template to reach a final reaction volume of 20 µl.

To minimize risk of contamination, always pipette templates in the following order: first NTC, then samples, last positive controls. Immediately close reaction vessels.

» **Note:** *Avoid creating bubbles in the final reaction mix and avoid touching the optical surface of the cap or sealing film without gloves. Both may interfere with fluorescence measurements. Centrifuge briefly if needed.* «

7.4. PCR Program

Program the real-time PCR instrument according to the manufacturer's instructions for allelic discrimination / genotyping experiments. Place the samples into the thermal cycler and run the following program:

Program			AB 7500 Fast, CFX96™, LightCycler® 480, and other Peltier heating-block based instruments	MIC, Rotor-Gene® 6000 (36-well & 72-well rotor)
Cycles	Temp	Time	Steps	Steps
1	95°C	3 min	Initial denaturation	Initial denaturation
40	95°C	15 sec	Denaturation	Denaturation
	60°C	1 min	Annealing/Extension –	Annealing/Extension –
			Data acquisition on FAM, HEX, ROX and Cy5 channels	Data acquisition on Green, Yellow, Orange and Red channels

8. Data Analysis / Interpretation of Results

The genotype of each sample is determined by calculating the ratio between signals recorded in the **HEX** or **Cy5 channel (normal)** and signals recorded in the **FAM** or **ROX channel (mutant)**. Most real-time PCR software automatically resolves data of two channels into clusters in a scatterplot. Data points plotted along the x- and y-axes correspond to normal and homozygous mutant genotypes, respectively. Data points clustered in the middle of the scatterplot represent heterozygous genotypes. The NTC appears in the lower left corner.

Controls / Samples	Amplification in channel				Genotype 677C>T / 1298A>C
	FAM Green	HEX Yellow	ROX Orange	Cy5 Red	
mpx WT-Control	NO	YES	NO	YES	normal 677C>T / normal 1298A>C
mpx HET-Control	YES	YES	YES	YES	heterozygous 677C>T / heterozygous 1298A>C
mpx MUT-Control	YES	NO	YES	NO	homozygous mutant 677C>T / homozygous mutant 1298A>C
NTC	NO	NO	NO	NO	----
Sample 1	YES	YES	NO	YES	heterozygous 677C>T / normal 1298A>C
Sample 2	YES	NO	NO	YES	homozygous mutant 677C>T / normal 1298A>C
Sample 3	NO	YES	YES	YES	normal 677C>T / heterozygous 1298A>C
Sample 4	NO	YES	YES	NO	normal 677C>T / homozygous mutant 1298A>C

Some instrument software needs manual threshold settings for accurate genotype calling.

Recommendations for Threshold Settings (C_q):

Set threshold value for the FAM and ROX channels just above the background fluorescent signal generated by the WT-Control (HEX-/Cy5-positive). Vice versa, set threshold value for the HEX and Cy5 channels just above the background fluorescent signal of the MUT-Control (FAM-/ROX-positive).

Samples crossing the threshold line beyond C_q 37 give invalid results and must be repeated.

To analyze acquired data, please follow your instrument software instructions.

9. Warnings and Precautions

- For *in vitro* diagnostics use only.
- Always use disposable powder-free gloves and wear suitable lab coat when handling specimens and reagents.
- Perform reaction setup in an area separate from nucleic acid preparation and PCR product analysis.
- Use pipettes dedicated for PCR setup only, use aerosol-guarded pipette tips.
- Use instrument-compatible reaction vessels with optically clear caps or sealers.
- Do not mix reagents from different lots.
- Do not use expired kits or kit components.

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