HLA-B1502 RealFast™ Assay

REF 7-630 / 7-633 Σ 100 / 32 reactions
-30°C -15°C CE IVD

ViennaLab Diagnostics GmbH
Gaudenzdorfer Guertel 43-45
A-1120 Vienna, Austria
Phone: (+43-1) 8120156-0
info@viennalab.com
www.viennalab.com

1. Intended Use

The HLA-B1502 RealFast™ Assay is a fast and accurate real-time PCR test for detection of the HLA-B1502 allele, a specific variant of the *human leukocyte antigen B (HLA-B)* gene that is strongly associated with carbamazepine-induced severe cutaneous adverse reactions. The kit is designed for genetic risk stratification of patients prior to initiation of carbamazepine therapy. HLA-B1502 positive patients must be excluded from carbamazepine treatment. The qualitative assay discriminates the presence or absence of HLA-B1502 in a human genomic DNA extract.

Reference sequence: NG_023187.1

2. Introduction

Carbamazepine is an anticonvulsant and mood-stabilizing drug commonly prescribed for multiple indications such as epilepsy, bipolar disorder, trigeminal neuralgia and chronic pain. In approximately 5% to 10% of individuals this medication can cause hypersensitivity reactions, of which Stevens-Johnson syndrome (SJS) and toxic epidermal necrolysis (TEN) are among the most severe complications. Although rare, SJS/TEN is responsible for a mortality rate as high as 50%, with fever, malaise and rapidly developing mucocutaneous blistering reactions as clinical manifestations. More extensive skin detachment and higher death rates have been observed in TEN. Numerous studies describe the strong association of the HLA-B1502 allele with carbamazepine-induced SJS/TEN in many Asian populations, where the prevalence of HLA-B1502 can be up to 15% or more in parts of China, Thailand, Malaysia, Indonesia, the Philippines and Taiwan. South Asians have an intermediate prevalence of HLA-B1502, whereas in Japan and Korea the allele is rare (<1%). Prospective HLA-B1502 screening of patients from at-risk populations therefore significantly reduces the incidence of hypersensitivity reactions due to carbamazepine.

3. Kit Contents

100 / 32 Rxn

RealFast™ 2x Genotyping Mix	1 vial	□ white cap	1000 / 320 µl
HLA-B1502 Assay Mix	1 vial	■ purple cap	550 / 550 µl
HLA-B1502 Positive Control	1 vial	■ green cap	75 / 75 µl
HLA-B1502 Negative Control	1 vial	■ red cap	75 / 75 µl

The RealFast™ 2x Genotyping Mix comprises HotStart Taq DNA polymerase and dNTPs in an optimized buffer system. The HLA-B1502 Assay Mix consists of gene-specific primers and dual-labeled hydrolysis probes for *HLA-B1502* and a control gene.

A positive and a negative control for HLA-B1502 are supplied with the kit.

The kit contains reagents for 100 / 32 reactions in a final volume of 20 µl each.

4. Storage and Stability

HLA-B1502 RealFast™ Assay is shipped on cooling blocks. On arrival, store the kit at -30 to -15°C. Alternatively, store at 2 to 8°C for short-term use within one month. The kit withstands up to 20 freeze/thaw cycles with no loss of activity. Avoid prolonged exposure to intense light. If stored correctly, the kit will retain full activity until the expiration date indicated on the label.

5. Product Description

5.1. Principle of the Test

The test is based on the fluorogenic 5' nuclease assay, also known as TaqMan® assay. Each reaction contains gene-specific primer pairs which amplify a 137 bp fragment of the HLA-B1502 gene and a 119 bp fragment of a control gene, the latter serving as PCR control. Further components are two dual-labeled, gene-specific hydrolysis probes which hybridize to the target sequence of the corresponding fragment. The proximity of the 5'-fluorescent reporter and 3'-quencher dye on intact probes prevents the reporter from fluorescing. During the extension phase of PCR the 5' – 3' exonuclease activity of the Taq DNA polymerase cleaves the 5'-fluorescent reporter from the hybridized probe. The physical separation of the fluorophore from the quencher dye generates a fluorescent signal in real-time, which is proportional to the accumulated PCR product.

In samples positive for HLA-B1502 both, the **FAM-labeled HLA-B1502** probe as well as the **HEX-labeled PCR control** probe bind to the appropriate gene fragment. A strong fluorescence signal is detected in the FAM channel (520nm) and in the HEX channel (556nm). In samples negative for HLA-B1502 only the HEX-labeled PCR control probe hybridizes to the complementary strand of the control gene fragment. A strong fluorescence signal is detected in the HEX channel and no or only a baseline signal in the FAM channel.

5.2. Real-time PCR Instrument Compatibility

The HLA-B1502 RealFast™ Assay is validated for use with the AB 7500 Fast instrument.

The kit is compatible with various common real-time PCR instruments capable of recording FAM and HEX fluorescence:

- ✓ AB 7500 Fast (Applied Biosystems®)
- ✓ AB StepOne™ (Applied Biosystems®)
- ✓ CFX96™ (Bio-Rad)
- ✓ LightCycler® 480 (Roche)
- ✓ MIC qPCR Cycler (bms)
- ✓ Rotor-Gene® 6000 (Qiagen)

» **Note:** RealFast™ Variant Detection QuickGuides for setting up and analyzing experiments on different types of instruments can be downloaded from www.viennalab.com.

When using AB StepOne™, set passive reference dye to "ROX!"

The kit is supplied **without ROX**. For use with real-time PCR instruments requiring high ROX for normalization of data (e.g. Applied Biosystems® instruments StepOne™, 7300, 7900/7900HT), add ROX at a final concentration of 1 µM to the 2x Genotyping Mix.

5.3. Assay Performance Specifications

Determination of **sensitivity** was performed on 58 alleles testing positive for the HLA-B1502 allele with Sanger sequencing and a CE-marked reference kit. The HLA-B1502 RealFast™ Assay determined all 58 alleles as positive, which equaled a true positive rate of 100%.

Determination of **specificity** was performed on 45 alleles testing negative for the HLA-B1502 allele with Sanger sequencing and a CE-marked reference kit. The HLA-B1502 RealFast™ Assay determined all 45 alleles as negative, which equaled a true negative rate of 100%.

Limit of detection: 2.5 ng genomic DNA (per reaction).

Recommended DNA concentration: 2 to 20 ng/µl genomic DNA.

6. Materials Required but not Supplied

Real-time PCR instrument with FAM (520 nm) and HEX (556 nm) filters, instrument-compatible reaction vessels, disposable powder-free gloves, vortexer, mini-centrifuge for 2.0 ml tubes, tube racks, set of calibrated micropipettes (0.5 – 1000 µl), sterile tips with aerosol-barrier filter, molecular grade water, DNA extraction system, freezer, biohazard waste container.

7. Experimental Protocol

7.1. DNA Extraction

DNA extraction reagents are **not supplied** with the kit.

DNA isolated from various specimens (e.g. whole peripheral blood, dried blood spots, buccal swabs or saliva) can be used. Ensure extracted DNA is suitable for amplification in terms of concentration, purity and integrity.

For accurate genotype calling, the DNA amount per reaction should be within the range of 10 to 100 ng for all samples.

7.2. PCR Controls

Always include a **No Template Control (NTC)** in each experiment to confirm absence of potential contamination. It is advisable to run the NTC (use PCR-grade water instead of DNA) in duplicate.

Always include the HLA-B1502 **Positive Control** as positive reference signal for your unknown samples and the HLA-B1502 **Negative Control** as negative reference signal for threshold setting in the FAM channel.

» **Note:** The Controls are potential sources of contamination. Make sure to handle them carefully. «

7.3. Preparation of HLA-B1502 RealFast™ Master Mix:

Gently vortex and briefly centrifuge all solutions after thawing. Set up PCR at room temperature. Prepare sufficient **Master Mix** for all your reactions (N samples + positive control + negative controls) plus at least one additional reaction to compensate for pipetting inaccuracies:

Component	per reaction	e.g. 24+1 reactions
RealFast™ 2x Genotyping Mix	10 µl	250 µl
HLA-B1502 Assay Mix	5 µl	125 µl
Master Mix	15 µl	375 µl

Dispense **15 µl Master Mix** into each well. Add **5 µl purified DNA or Control** template to reach a final reaction volume of 20 µl.

To minimize risk of contamination, always pipette templates in the following order: first NTC, then samples, last controls. Immediately close reaction vessels.

» **Note:** Avoid creating bubbles in the final reaction mix and avoid touching the optical surface of the cap or sealing film without gloves. Both may interfere with fluorescence measurements. Centrifuge briefly if needed. «

7.4. PCR Program

Program the real-time PCR instrument according to the manufacturer's instructions for quantitation experiments with two targets / reporter dyes. Place the samples into the thermal cycler and run the following program:

AB 7500 Fast, StepOne™, CFX96™, LightCycler® 480
and **other Peltier heating block-based instruments:**

Cycles	Temp	Time	Steps
1	95°C	3 min	Initial denaturation
40	95°C	15 sec	Denaturation
	60°C	1 min	Annealing/Extension – Data acquisition on FAM and HEX channel

MIC qPCR Cycler, Rotor-Gene® 6000*):

Cycles	Temp	Time	Steps
1	95°C	3 min	Initial denaturation
40	95°C	15 sec	Denaturation
	60°C <i>*)for 36-well rotor: 56°C</i>	1 min	Annealing/Extension – Data acquisition on Green and Yellow channel

8. Data Analysis / Interpretation of Results

The presence or absence of the HLA-B1502 allele is defined by whether there is a signal in the **FAM channel** or not. Successful PCR can be verified by amplification of the control gene detected in the **HEX channel** (PCR control). Thus, genomic DNA samples positive for HLA-B1502 as well as the HLA-B1502 Positive Control show amplification in both, the HEX and FAM channel. HLA-B1502 negative samples and the HLA-B1502 Negative Control show amplification in the HEX channel only. Fluorescent levels and corresponding amplification curves are automatically displayed in amplification plots in the real-time PCR software.

Sample Type	Amplification in FAM channel (520 nm)	Amplification in HEX channel (556 nm)
HLA-B1502 positive	YES	YES
HLA-B1502 negative	NO	YES
HLA-B1502 Positive Control	YES	YES
HLA-B1502 Negative Control	NO	YES
NTC	NO	NO


Some instrument software may need manual threshold settings for accurate analysis.

Recommendations for Threshold Settings (C_q):

Set threshold value for the FAM channel just above the background fluorescent signal generated by the HLA-B1502 Negative Control.

Samples crossing the threshold line beyond C_q 37 give invalid results and must be repeated.

To analyze acquired data, please follow your instrument software instructions.

9. Warnings and Precautions

- For *in vitro* diagnostics use only.
- Always use disposable powder-free gloves and wear suitable lab coat when handling specimens and reagents.
- Perform reaction setup in an area separate from nucleic acid preparation and PCR product analysis.
- Use pipettes dedicated for PCR setup only, use aerosol-guarded pipette tips.
- Use instrument-compatible reaction vessels with optically clear caps or sealers.
- Do not mix reagents from different lots.
- Do not use expired kits or kit components.

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