GEN^XTRACT™ Blood DNA Extraction System

REF

2-014

IVD

100 extractions

2-8°C

CE

1. Lysis Solution	4x 50 ml
2. GEN ^X TRACT™ Resin	4x 5 ml
3. Instructions For Use	1

ViennaLab Diagnostics GmbH

Gaudenzdorfer Guertel 43-45
A-1120 Vienna, Austria
Phone: (+43-1) 8120156-0
Fax: (+43-1) 8120156-19
info@viennalab.com

www.viennalab.com

Instructions for use

I. INTENDED USE

GEN^XTRACT™ Blood DNA Extraction System is designed to be used for rapid and easy extraction of DNA from blood as template for *in vitro* amplification (PCR). *For human in vitro diagnostics.*

Store all reagents at 2-8°C when not in use !

II. MATERIALS REQUIRED BUT NOT SUPPLIED

- Microtubes (1.5 ml with screw cap)
- Adjustable micropipettes
- Adjustable microcentrifuge capable of 3,000-12,000 rpm (1,000-12,000 x g)
- Incubator (e.g. heating block, water bath) capable of 56°C and 98°C (± 2°C)

III. ASSAY PROCEDURE

1. DNA Isolation from Whole Blood

Use fresh or frozen blood with EDTA or citrate anticoagulant; avoid blood containing heparin. Do not store blood for more than 3 days at ambient temperature or more than 1 week at 2-8°C before use. Blood which has been kept frozen for more than one year, or gone through more than three freeze-thaw cycles is unsuitable to be used in this procedure.

Bring blood samples to room temperature. Mix well by carefully inverting blood collection tubes several times. Repeat mixing each time before withdrawing an aliquot of blood.

Allow Lysis Solution and GEN^XTRACT™ Resin to reach room temperature.

- Pipette **100 µl blood sample** into a 1.5 ml microtube with screw cap.
- Add **1 ml Lysis Solution**, close tube and mix by inverting several times.
- Let stand for **15 min.** at room temperature.
- Centrifuge for **5 min.** at **3,000 rpm** (approx. 1,000 x g) in a microcentrifuge.
- Remove and discard the upper (top) 1 ml of supernatant.
- Add **1 ml Lysis Solution**, close tube and mix by inverting several times.
- Centrifuge for **5 min.** at **12,000 rpm** (approx. 12,000 x g) in a microcentrifuge.
- Remove and discard the supernatant except for approx. 50 µl of a visible, soft pellet.
- Resuspend GEN^XTRACT™ Resin by swirling the bottle thoroughly.
- Add **200 µl GEN^XTRACT™ Resin** to the pellet. Close tube and vortex for 10 sec.
⚠ GEN^XTRACT™ Resin sediments quickly. Repeat resuspension each time immediately before removing another aliquot.
- Incubate for **20 min.** at **56°C**. Vortex for 10 sec.
- Incubate for **10 min.** at **98°C**. Vortex for 10 sec.
- Centrifuge for **5 min.** at **12,000 rpm** in a microcentrifuge. Cool on ice.

The resulting supernatant contains DNA template suitable for immediate use in PCR. For further storage, the supernatant should be transferred into a fresh tube and kept refrigerated (2-8°C; up to one week) or frozen at -20°C.

2. DNA Isolation from Dried Blood Spots

Collect blood drops (from a finger or heel prick, or EDTA anticoagulated blood) onto a ring marked on a Whatman 903[®] sample collection card. Allow the blood to air-dry.

For long term archiving, store blood cards at 2-8°C in a resealable bag with a desiccant pouch. Allow Lysis Solution and GEN^XTRACT™ Resin to reach room temperature.

- Place **two 3 mm punches** into a 1.5 ml microtube with screw cap.
- Add **1 ml Lysis Solution**, close tube and mix by inverting several times.
- Let stand for **10 min.** at room temperature and mix by inverting several times.
- Centrifuge for **1 min.** at **12,000 rpm** (approx. 12,000 x g) in a microcentrifuge.
- Remove and discard the supernatant completely.
- Add **1 ml Lysis Solution**, close tube and mix by inverting several times
- Let stand for **10 min.** at room temperature and mix by inverting several times.
- Centrifuge for **1 min.** at **12,000 rpm** in a microcentrifuge.
- Remove and discard the supernatant completely.
- Resuspend GEN^XTRACT™ Resin by swirling the bottle thoroughly.
- Add **200 µl GEN^XTRACT™ Resin** to the punches. Close tube and mix gently by tapping the bottom of the tube. Punches should be completely submerged by the resin.
⚠ GEN^XTRACT™ Resin sediments quickly. Repeat resuspension each time immediately before removing another aliquot.
- Incubate for **20 min.** at **56°C**.
- Incubate for **10 min.** at **98°C**.
- Centrifuge for **5 min.** at **12,000 rpm** in a microcentrifuge. Cool on ice.

The resulting supernatant contains DNA template suitable for immediate use in PCR. For further storage, the supernatant should be transferred into a fresh tube and kept refrigerated (2-8°C; up to one week) or frozen at -20°C.

IV. QUALITY CONSIDERATIONS

- A thorough understanding of the procedure outlined here, and precise laboratory equipment and techniques are required to obtain reliable results. Use of GEN^XTRACT™ Blood DNA Extraction System for human *in vitro* diagnostics needs to be limited to appropriately trained personnel.
- Do not use GEN^XTRACT™ Blood DNA Extraction System beyond the expiration date printed on the outside of the kit box. Do not mix reagents from different lots.
- Avoid microbial contamination and cross-contamination of reagents or samples by using sterile disposable pipette tips throughout. Do not interchange bottle caps.

V. SAFETY

- Do not drink, eat, smoke, or apply cosmetics in designated work areas. Wear laboratory coats and disposable gloves when handling specimens and kit reagents. Wash hands thoroughly afterwards.
- Handle specimens as if capable of transmitting infectious agents. Thoroughly clean and disinfect all materials and surfaces that have been in contact with specimens. Discard all waste associated with clinical specimens in a biohazard waste container.
- Adhere to all local and federal safety and environmental regulations which may apply.

VI. TROUBLESHOOTING

Advice on troubleshooting may be obtained by contacting ViennaLab through the local distributor or directly at techhelp@viennalab.com.

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