


StripAssay[®] Evaluator

StripAssay[®] Evaluator Software

Our Evaluator software permits one-step, automated analysis of up to 20 teststrips. The software converts the color intensity of the StripAssay[®] bands into digital data that can be stored electronically. The software allows precise and effortless analysis of complex StripAssays[®].



Easy Handling Test type is automatically recognized via barcodes

Easy Interpretation Results are automatically analyzed and displayed in comprehensive form

Easy Safe-Guard Unusual results trigger an inbuilt warning system

Easy Data Storage Test results can be archived, exported and cross-referenced